

Reanibex 200

Automated External Defibrillator

2 versions:

- 1. Icons are displayed
- 2. Icons or ECG are displayed

Biphasic up to 200 J

Universal AED algorithm for adult and paediatric patients

Lithium battery with more than 400 discharges of capacity

Data transmission with IrDA

Automatic daily selftest

CPR Metronome with feedback

Reanibex 200.

Portable & lightweight



FEATURES

- Biphasic technology.
- Universal AED algorithm for adult and paediatric patients.
- Energy delivered: 200 Joules for adult and 50 Joules for paediatric.
- Audible and visual prompts to guide users through complete operation sequence.
- Configurable CPR protocol.
- Lithium battery with more than 400 discharges of capacity.
- Configured according to ERC/ AHA guidelines 2015 and its review of 2017.
- Automatic daily selftest ensures the device continual readiness.
- Portable and lightweight for easy transport in any situation.
- Rugged and durable design for any environment.
- Optional PC software, to store and review ECG data.


ACCESORIES

- Adult and paediatric disposable electrodes for defibrillation.
- Carrying case (1).
- Ambulance bracket (EN1789:2007 certified) (2).
- Non rechargeable battery (3).
- "REANIBEX Data Manager" Software application to manage and organize all the data gathered (4).
- Wallcase with optional alarm system. Other models available with Communication systems (GSM, GPRS), CCTV camera...).











Bexen cardio